

## Capillary HPLC:

### Fast PLRP-S Separation of A Small Synthetic Peptide




App No 112

#### Small Synthetic Peptide

**Column:** 3 $\mu$  PLRP-S 100Å  
**Dimensions:** 0.3 x 50 mm  
**Order No.:** S-161-0305  
**Mobile Phase:** A: 0.1% TFA in MeCN-water (1:99)  
 B: 0.1% TFA in MeCN-water (50:50)  
**Gradient:** 40-80% in 5 mins  
**Flow Rate:** 5  $\mu$ L/min  
**Temp.:** 30C  
**Detector:** UV @ 220 nm

Ac-Lys-Tyr-Ala-Leu-Lys-Ala-  
 Leu-Lys-Gly-Leu-Lys-acid



#### Comments:

Small peptides usually need the high surface area of  $\leq 100\text{\AA}$  particles to obtain retention.

Trifluoroacetic acid (TFA) is usually the ion-pairing agent of choice when separating peptides by HPLC-UV.

#### ORDERING INFORMATION: HotSep® PLRP-S, 3 $\mu$ , 100Å

ID	3 cm	5 cm	10 cm	15 cm	Guard	Guard/5pk
1.0 mm	S-161-1003	S-161-1005	S-161-1010	S-161-1015	G-161-10-1	G-161-10-5
0.5 mm	S-161-0503	S-161-0505	S-161-0510	S-161-0515	G-161-05-1	G-161-05-5
0.3 mm	S-161-0303	S-161-0305	S-161-0310	S-161-0315	G-161-03-1	G-161-03-5
0.1 mm	---	---	S-161-0110	S-161-0115	---	---
75 $\mu$ m	---	---	S-161-00710	S-161-00715	---	---