



## Microbore HPLC: Separation of Fat Soluble Vitamins



App No 127 **Fat Soluble Vitamins**

**Column:** 3 $\mu$  ProntoSIL C18 100Å  
**Dimensions:** 0.5 x 150 mm  
**Order No.:** S-132-0515  
**Eluent:** MeCN  
**Flow Rate:** 10  $\mu$ L/min  
**Temp.:** 30C  
**Detector:** UV @ 280 nm



### ORDERING INFORMATION: HotSep® ProntoSIL C18, 3 $\mu$ , 100Å

ID	3 cm	5 cm	10 cm	15 cm	Guard	Guard/5pk
1.0 mm	S-132-1003	S-132-1005	S-132-1010	S-132-1015	G-132-10-1	G-132-10-5
0.5 mm	S-132-0503	S-132-0505	S-132-0510	S-132-0515	G-132-05-1	G-132-05-5
0.3 mm	S-132-0303	S-132-0305	S-132-0310	S-132-0315	G-132-03-1	G-132-03-5
0.1 mm	---	---	S-132-0110	S-132-0115	---	---
75 $\mu$ m	---	---	S-132-00710	S-132-00715	---	---